MACROMAN



Stable, safe & easy to pick up

- Integrated bench stand adjustable for left and right handed users.
- > **Keep your MACROMAN handy** and ready for use, without risk of contamination.

Soft & fatigue free pipetting

- > **Soft-touch material** for maximum finger comfort.
- > Large dosing button to decrease pushing efforts.





MACROMAN

Optimal comfort

Integrated, adjustable stand for easy handling when pipette

Large & soft-touch material dosing button

Light weight

Quiet operation & vibration free

Ergonomic-shape handle with grip



All liquid contact parts are autoclavable

Easy-to-replace filter

Replaceable batteries

Performance

For all plastic & glass pipettes from 1 to 100mL

Continuously adjustable speed: 25mL in 4sec

Environmentally friendly NiMH battery

for **8 hours operation**

Precise & smooth dosing button control: drop-by-drop, gravity mode and motorized speed.

Convenient accessories

Horizontal or vertical & magnetic charging stand (with complete mounting kit) for busy laminar flow cabinet.

Unique adaptor for D5000 and D10mL DIAMOND® tips.



MACROMAN

Weight	194g
Battery	2x 1.2V AAA NiMH, removable, 8-hour operation, rechargeable during use, discharged & charging light.
Autoclavability	Nose piece, Silicone pipette mount & filter rubber.
Filter	Replaceable, 0.45μm or 0.20μm standard & sterile.
Pipette types	Glass or plastic from 1mL to 100mL
Power supply	100-240V / 50-60Hz/ EU, US, UK, JAP or AUS plug version
Certification	CE, EMC, RoHS, WEEE
Contents	Complete MACROMAN + holder +1 extra 0.45 µm sterile filter + identification labels + EU, US, UK, JAP or AUS plug version

Description	Reference Number
MACROMAN, EU plug	F110751
MACROMAN, US plug	F110752
MACROMAN, UK plug	F110753
MACROMAN, JAP plug	F110754
MACROMAN, AUS plug	F110755
Charging stand	F1077501
0.45µm standard filter	F1077504
0.45µm sterile filter	F1077525
0.20μm standard filter	F1077503
0.20μm sterile filter	F1077526
Nose piece	F1077506
Filter rubber	F1077502
Pipette mount	F1077505
Adapter for D5000 and D10mL DIAMOND® Tips	F161287

Distributed by

Gilson provides 1-year warranty after the date of purchase

www.gilson.com | sales@gilson.com | service@gilson.com | training@gilson.com

Gilson, Inc. World Headquarters | 3000 Parmenter Street | P.O. Box 620027 | Middleton, WI 53562-0027 USA Tel: (1) 800-445-7661 or (1) 608-836-1551 | Fax: (1) 608-831-4451

Gilson S.A.S. | 19, avenue des Entrepreneurs | BP 145, F-95400 Villiers-le-Bel, France Tel: (33-1) 34 29 50 00 | Fax: (33-1) 34 29 50 20

LT800498E - © 2009 Gilson S.A.S. - All rights reserved - Printed in France, June 2009. Gilson, Inc. and Gilson S.A.S are ISO 9001 certified. Specifications subject to change without notice.

